Technical Description Sheet

HYGIEN'AIR

Hygien'Air is a nebulizing spray machine that rapidly disinfect defined spaces and surfaces by using high technical solutions



Diffusions 1ml/m³ • 5ml/m³ • 7ml/m³



Contact us: 00972-9-885-8896 +972-54-482-0549 00972-9-955-9309 WWW.ouriel-medical.com

Ouriel Medical Equipment for Hospitals and Research Laboratories

HYGIEN'AIR

Technical Information

Hygien'Air is a light, ergonomic nebulizing device Weight: 4,5 kg It disinfects volumes from 5m³ to 165m³ Voltage: 230V or 115V Diffusion time: 13,2 sec./m3

Existing Devices

- ✓ Hygien'Air 1M: diffusion of 1ml/m3 for insecticide The choice can depend on the applications and sectors
- ✓ Hygien'Air 5M: diffusion of 5ml/m3
- ✓ Hygien'Air 7M: diffusion of 7ml/m3

Applications

Hygien'Air can be used with appropriate disinfectants solutions, insecticides or desodorisant Hygien'Air can be used in the following areas where hygiene is important: -Medical and para-medical fields - Ambulances -Pharmaceutical industry - Laboratories and clean-rooms -Veterinary field - Military field - Agro-food industry - Schools...

Advantages

- Efficient : decontaminate and disinfect spaces and areas that are difficult to reach manually
- Homogeneous disinfection action
- ✓ Safe for the users: automated device that switches on automatically after 90 sec. or 600 sec.
- Simple and fast: Easy to program the device with the display and to set it according to the desired quantity of the solution to be used
- Dry nebulization, no humidity on the surfaces > no risk of corrosion
- ✓ Ecological disinfection, environment-friendly
- ✓ Cost-effective application



Contact us: ☎ 00972-9-885-8896 ■ +972-54-482-0549 ■ 00972-9-955-9309 ⊠ info@ouriel-medical.net WWW.ouriel-medical.com